**HAYER & BOECKER**

**DIE DRAHTWEBER**

## **HAYER RPT**

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## *Index*

		<i>Page</i>
<b>1.</b>	<b><i>Introduction</i></b>	<b>13</b>
<b>2.</b>	<b><i>Transport</i></b>	<b>13</b>
<b>3.</b>	<b><i>Proper use</i></b>	<b>14</b>
<b>4.</b>	<b><i>Principles of operation</i></b>	<b>14</b>
<b>5.</b>	<b><i>Scope of delivery</i></b>	<b>15</b>
<b>6.</b>	<b><i>Technical data</i></b>	<b>16</b>
<b>7.</b>	<b><i>Operational safety</i></b>	<b>16</b>
<b>8.</b>	<b><i>Setting up the device</i></b>	<b>17</b>
<b>9.</b>	<b><i>Sample division</i></b>	<b>18</b>
<b>10.</b>	<b><i>Control</i></b>	<b>20</b>
<b>11.</b>	<b><i>Guarantee</i></b>	<b>21</b>
<b>12.</b>	<b><i>Maintenance and repairs</i></b>	<b>21</b>
	<b><i>Appendix: EC-Attestation of conformity, Table of available test sieves and accessories</i></b>	

## 1. Introduction

**Do not set up the machine or put it into operation until after you have read the following Operating Manual thoroughly.**

The Operating Manual contains important information on safe, correct, and efficient machine operation. Using the machine in accordance with these instructions will help to avoid any dangers, minimize repairs and downtime, improve machine reliability, and to prolong the service life of the machine.

When operating this machine, it is essential that you observe the general safety instructions and the regulations for the prevention of accidents as well as the protective measures recommended by the professional and trade associations, which may vary depending on country and/or professional association.

**Adhere to all safety instructions contained in this operating manual.**

**The operating instructions must always be available at the place where the machine is in use.**

When operated according to its intended purpose, the design of the machine guarantees that it can be operated and maintained without risk of injury.

**Any damage resulting from non-compliance with the instructions in the Operating Manual is not covered by the guarantee.**

The purpose of all figures and drawings in this Operating Manual is to provide a general overview of the machine. They are not intended to show details of the design.

The technical specifications and dimensions are without obligation. No claims can be derived based on these specifications.

We reserve the right to make technical changes resulting from further development of the machine covered by this operating manual without changing the manual.

All rights to this operating manual, to the attached drawings and other documentation as well as the power of disposal such as the right to copy and transmit the documents lie with Haver & Boecker, including applications for protective rights.

### Symbols used in this operating manual

**Caution!**  
**Warning, hazardous area**  
**Adhere to instructions in operating manual**

**Caution!**  
**Supply voltage**

**Caution!**  
**Flammable substances**

**Wear eye protectors**

## 2. Transport

With the exception of the divider core and the hopper, the HAVER RPT is delivered fully assembled. Unless otherwise agreed, the machine is packed in accordance with the HPE packaging guidelines laid down by the "Bundesverband Holzmittel, Paletten, Exportverpackungen e.V." (the national association of wooden instruments, pallets, export packaging) and the "Verein Deutscher Maschinenbauanstalten" (the society of German engineering shops)

Observe the relevant national and international safety prevention regulations currently in effect. Damage due to incorrect transportation does not qualify for replacement or entitlement under the guarantee.

**Carry the machine by its base plate and not by the feeder!**

### 3. Proper Use

The HAVER RPT is a rotating sample reducer equipped with a vibration feeder.

 The machine is constructed in accordance with the best available technology and recognized safety-related regulations. Nevertheless, it is possible that operating the machine may endanger the user or other persons, or have a detrimental effect on the machine or other property.

The machine may be used only if in perfect working condition and in the proper manner with attention to safety and risks in accordance with these operating instructions. In particular, any faults that may affect safe operation must be corrected immediately!

**The rotating sample reducer** is used for a representative sample division in a laboratory. It divides dry bulk materials or particle suspensions of a feed rate in similar quantities with each partial quantity being representative for the total sample.

Inorganic or organic samples are prepared for analysis, quality control or materials testing in a way that after analysing a single sample, statements on physical or chemical characteristics of the total sample will be reliable.

**The vibration feeder** is used for continuous or discontinuous feeding of the rotating sample reducer. It uniformly feeds fine-grained to coarse free-flowing solids.

Materials exposed to the flow of bulk materials are made of stainless steel.

The feed rate may be set from 1g/min to 2.5kg/min.

The maximum feeding size of a single particle is an edge length of approx. 10 mm.

### 4. Principle of Operation

The HAVER RPT combines two different principles of division in one machine in a way that the advantages of both may be optimally used:

The bulk material flow led through a hopper impacts against a cone. The first division into single flows takes place while sliding its cylindrical surface. This method is known as "coning" is successfully used above all for dividing larger quantities.

The HAVER RPT directs the expanded bulk material flow in single channels at the lower end of the cone's cylindrical surface where it is then collected in a laboratory glass jar. Turning the cone with a relative high rotational speed avoids inhomogeneity in the bulk material flow or segregations during feeding. Rotating the cone equalises the chance of all amounts of bulk material falling into a particular glass jar. Additionally, the separated bulk material flow is accelerated outwards by rotation and flows more easily through the guide channels into the collecting glass jar. The guide channels permeate each other. The sharp dividing lines between them partition the bulk material flow and minimise deposits. The high rotational speed of the circumferential dividing cone causes a very high dividing ratio, i.e. each sample is composed of a very high number of individual samples. A large dividing ratio, however, is one of the most important prerequisites for a good and representative sample dividing.

The stainless steel chute of the feeder is vibrated by an electromagnet. A height-adjustable stainless steel hopper plunges into this feeder fixed to a supporting column.

The materials conveyed are filled into the hopper. When the feeder is switched on, the clearance between hopper and feeder defines the amount of outflowing bulk material.

The electronic control defines the vibration amplitude of the vibrating feeder and thus the amount and flow rate of the material being fed.


The housing contains an electronically controlled electromagnet. Four permanently elastic springs support the anchor of this magnet with a work plate attached to it so that when the magnet is activated the tension of the work plate is increased enabling it to spring back when switched off. Housing and magnet on one side together with anchor and work plate on the other comprise a vibrating system. Its natural frequency, for example, is influenced by the mass of the bulk material.

The desired feed rate is reached under all operating conditions; processor controlled electronics provide a reproducible vibration amplitude by tuning the frequency of the frequency generator to the natural frequency of the vibration feeder.

The interface integrated into the feeder control allows the feeder flow to be stopped by an external signal, so that for example the feed rate may be controlled by a connected direct current (range of 5-30V).

"Remote control" of the vibration feeder facilitates the process operation or routine work in the laboratory significantly while at the same time ensuring reproducible process steps.


## 5. Scope of Delivery


HAVER RPT


Feeder hopper

Feeder

One divider core including laboratory glass jars


Divider core cover

Hopper  $\varnothing$  10 mm for sample division at the rate of 1 : 8 and 1 : 10;  
 $\varnothing$  5 mm for sample division at the rate of 1 : 30

Supplementary hoppers are available for an optimum adaption to the respective particle size.

1 connection cable, 1.5 m long, with EURO plug and device connector

## 6. Technical Data of the HAVER RPT

Operating voltage:	230 V or 110 V with transformer
Frequency:	50 Hz / 60 Hz
Weight:	30 kg
Dimensions:	660 x 620 x 385 mm (Height x Width x Depth)
Feed rate per minute:	1 g/min to 2.5 kg/min
Materials exposed to the bulk material flow:	Food-safe aluminium, stainless steel, polyoxymethylene, glass

Divider head with division ratio	Maximum quantity for one run	Maximum permissible particle size
1:8	4000 ml	10 mm
1:10	2500 ml	10 mm
1:30	300 ml	2.5 mm

## 7. Operational safety

### General safety information

- Carefully read the operating manual.
- The unit may be operated only for its intended areas of application.
- Always use original accessories and spare parts. Failure to comply may jeopardize machine protection.
- Do not use damaged accessories.
- Device operators must be familiar with the contents of the operating manual. The operating manual must always be available at the place where the machine is being operated.
- Do not render the safety equipment inoperable.
- Unauthorised modifications made to the machine will void the declared conformity with European directives and void the warranty.
- Strictly adhere to rules of safe working conduct and procedures.

### Operating personnel

- The unit may only be operated by authorized persons and maintained and repaired only by trained personnel.
- The rotating sample reducer is not to be operated by persons suffering from physical disorders, under the influence of drugs, alcohol, or extreme fatigue.

### **Special safety information**

- "Out of service" ⇒ *turn off the main switch!*
- *Start running the sample reducer only with glass jars or collecting jars securely attached!*
- *Open the rotating sample reducer only after it has come to a full stop!*
- *Do not reach into the device while it is running!*

- *Always wear safety goggles.*
- *Set up the rotating sample reducer only in indoors.*
- *Ambient temperature range: 0 - 40°C.*

- *Bench-top installation.*
- *Danger! Electrical voltage!*
- *Slight risk of crushing from rotating glass jars: Attach cover!*
- *Exercise caution while using flammable or toxic substances!*

- *When using combustible or harmful samples or suspending liquids, always follow the safety regulations (MAK values) currently in effect and, if necessary, set up the rotating sample reducer in a ventilated safety zone.*

## **8. Setting up the device**

### **Setting up the Haver RPT**

- *Place the device on a solid, horizontal surface.*
- *The rubber feet provide a safe stand without the need of fastening the device to the working surface.*
- *Pay particular attention to ensuring easy access.*
- *The device may be operated only indoors.*
- *Ambient air may not contain any electrically conductive dust.*
- *Ensure that the room temperature is between 5 and 40°C.*
- *Operate at an altitude of up to 2000m*
- *Maximum relative humidity 80% for temperatures of up to 31°C, decreasing linear down to 50% relative humidity at 40°C*
- *Degree of soiling 2 according to IEC 664*

### **Attaching the divider core**

- *Attach the divider core to the drive axle protruding out of the housing.*
- *Check: The divider core should rotate smoothly by hand on the friction clutch.*

### **Electrical power supply**


- *Check for compliance of voltage and current on the nameplate with the values of the designated mains.*

## 9. Sample division

### 9.1. General preparations

#### Screw fitting of laboratory glass jars

- Align the mark on the glass jar with the mark on the glass jar holder
- Insert the glass jar into the cut-out of the holder
- Screw in the glass jar tightly clockwise  $\frac{1}{4}$  turn

#### Division ratio of 1:8

- Screw eight laboratory glass jars into the divider core (eight 250-ml jars or eight 500-ml jars depending on the quantity of the sample)
- Place the cover on the divider core
- Insert the hopper

#### Division ratio of 1:10


- Screw ten 250-ml laboratory jars into the divider core
- Place the cover on the divider core
- Insert the hopper
- 500-ml jars cannot to be used

#### Division ratio of 1:30

- Screw three laboratory glass jars into the divider core (three 15-ml jars, 20-ml jars, or 30-ml jars depending on the quantity of the sample)
- Insert and lock three catch trays (pull down the pin, insert the tray, release the pin)

#### Insert the feeder


- Insert the feeder into the guide so that the **run-out centred above the hopper**
- Adjust the amount of bulk material flowing out of the hopper to the feeder by altering its height at the supporting column. In case of dusty materials, do not select a clearance that is too large (approx. 3 mm) to avoid an excessive amount of dust and/or overflowing.

**Note: Do not fill with any material before switching on the sample divider!**


## 9.2. Dividing solids

- Insert hopper in position 1 (for solids)
- Switch on sample divider (wait for final speed of approx. 100 rpm)
- Feed by hand:
  - Fill the sample in the hopper evenly.
- Add material using the vibration feeder:
  - Turn feed rate (speed of material flow) to the left stop at the control potentiometer
  - Fill the sample into the feeder hopper.
  - Slowly turn the potentiometer clockwise until the desired feed rate is reached

## 9.3. Dividing suspensions

- Remove the feeder
- Insert hopper in position 2 (for suspensions)
- Switch on sample divider (wait for final speed of approx. 100 rpm)
- Pour the sample manually and evenly in the hopper.

## 9.4. Removing samples

- Turn off the device at the end of the dividing process
- Screw off the glass jars and remove the divided samples

## 9.5. Cleaning

- If necessary, remove divider core from the drive axle for cleaning (sliding disk is glued to the divider core)
- Keep dust off the drive axle of the gear motor and oil it lightly with machine oil once a month.
- Attach the divider without using force

### Normal soiling

- Vacuum out the hopper, divider core, and sample glass jars

### Heavy soiling

- Rinse hopper, divider, and sample glass jars with a soft bottle brush under running water

### Drying

- Rinse and air-dry hopper, divider, and sample jars with alcohol (ethyl alcohol) after wet cleaning


## 9.6. Maintenance

- Clean sample divider regularly
- Lightly oil the drive axle with machine oil once a month
- The ball bearings are lubricated for life
- No further maintenance is required

## 10. Control

*Rocker power switch:  
The switch will light when the device is ON  
The rotating sample reducer rotates*

*Control the feed rate in the feeder  
Infinitely variable vibration amplitude  
Setting the potentiometer 0 stops the feeder*

## 11. Guarantee

Each device is set and all functions are tested prior to delivery.

The guarantee period is **24 months**.

Our guarantee covers only defects based on manufacturing faults or material flaws. Guarantee does not apply to damage incurred during shipment. Guarantee does not include repairs caused by improper use or by operational wear. Guarantee is voided in the event the device is modified, improperly maintained, or otherwise altered.

## 12. Maintenance and repairs

**Please note that the HAVER-RPT may only be repaired, serviced, and inspected according to national regulations (e.g. according to VBG 4) by authorised electricians as improper repairs may lead to serious hazards to the user.**

**If a part of the machine and its associated equipment is replaced or modified, this part has to be retested according to EN 60 204-1.**

**We recommend that you contact our Service Dept. concerning any malfunctions of the HAVER rotating sample reducer.**

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**The HAVER rotating sample reducer should be inspected and serviced at our factory every two years at the latest.**

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